

HY2000 Transfection Reagent

保存：2-4°C保存一年。（避免冷冻）

产品说明

HY2000 是一种新型的阳离子脂质体转染试剂。适合于将核酸(DNA 和 RNA)转染入真核细胞，具有低细胞毒性；对多种类型的细胞和培养板都具有高转染效率；转染时血清的存在不影响转染效率的优点。

适用范围：贴壁细胞和悬浮细胞（哺乳动物细胞系）的转染。

质粒 DNA 的转染

对大多数细胞来说，DNA(μg)与 HY2000 (μl)的比例为 1:2~1:3。转染时高的细胞密度可以得到高的转染效率和表达水平，并能减少细胞毒性。

1. 以 24 孔板为例

贴壁细胞：转染前一天，用 500 μl 不含抗生素的培养基接种 $0.5\text{-}2 \times 10^5$ 细胞，使之第二天能达到 70-

90%汇合。悬浮细胞：在准备 DNA-HY2000 复合物之前，用 500 μl 不含抗生素的培养基接种 $4\text{-}8 \times 10^5$ 细胞即可。

2. 对每个转染样品，进行以下操作

a. 在 eppendorf 管里分别加入 50 μl Opti-MEM I ReLipced Serum Medium 和 0.8 μg DNA，轻柔混匀，制成 DNA 稀释液。

b. 在另一个 eppendorf 管里分别加入 50 μl Opti-MEM I ReLipced Serum Medium 和 2.0 μl HY2000（注意用前先混匀），轻柔混匀，制成 HY2000 稀释液，室温静置 5 分钟。

c. 将 DNA 稀释液和 HY2000 稀释液混合，轻柔混匀，室温静置 20 分钟，形成 DNA-HY2000 复合物。DNA-HY2000 复合物在室温下可稳定存在 6 小时。

3. 将 DNA-HY2000 复合物加入到接种好的细胞中，将培养板轻轻地前后摇动，使复合物分散均匀。

4. 在 37°C CO₂ 培养箱中培养 4-6 小时后更换培养基，继续培养 18-48 小时。

5. 如果要筛选稳定细胞株，则在转染 24 小时后将细胞按照 1:10 或更高的比例接种到新鲜培养基中，第二天加入选择性培养基进行筛选。

质粒 DNA 转染的优化 为达到最高的转染效率和降低细胞毒性的影响，可以对 DNA 和 HY2000 的比例以及细胞密度进行优化，一般在 1:0.5~1:5 的范围内优化 DNA (μg) 和 HY2000 (μl) 的比例。

不同细胞培养板中转染时培养基、核酸及 HY2000 用量

细胞培养板	每孔面积	培养基用量		DNA 转染		siRNA	
		铺板培养基用量	稀释培养基用量				
96-well	0.3 cm ²	100 μl	$2 \times 25 \mu\text{l}$	0.2 μg	0.5 μl	5 pmol	0.25 μl
24-well	2 cm ²	500 μl	$2 \times 50 \mu\text{l}$	0.8 μg	2.0 μl	20 pmol	1.0 μl
12-well	4 cm ²	1 ml	$2 \times 100 \mu\text{l}$	1.6 μg	4.0 μl	40 pmol	2.0 μl
6-well	10 cm ²	2 ml	$2 \times 250 \mu\text{l}$	4.0 μg	10 μl	100 pmol	5 μl
60-mm	20 cm ²	5 ml	$2 \times 0.5 \text{ ml}$	8.0 μg	20 μl	200 pmol	10 μl
10-cm	60 cm ²	15 ml	$2 \times 1.5 \text{ ml}$	24 μg	60 μl	600 pmol	30 μl

HY2000 Transfection Reagent

Description

HY2000 is a newly developed and proprietary reagent for the transfection of nucleic acids into eukaryotic cells.

HY2000 has the following advantages:

The highest transfection efficiency in many cell types and formats.

DNA-HY2000 complexes can be directly added to cells in culture medium (with or without serum).

It is not necessary to remove DNA-HY2000 complexes or change medium following transfection.

The complexes can be removed after 4-6 hours by replacing with refresh medium (optional)

Contents and Storage

HY2000 is supplied in liquid form at a concentration of 1mg/ml. Store at 4°C. **DO NOT FREEZE.**

Product Qualification

HY2000 has been extensively tested by transfection of HEK293 cells with an EGFP reporter containing plasmid. HY2000 is free of microbial contamination.

Important Guidelines

Follow these guidelines when performing transfections:

1. The ratio of DNA (in μg) : HY2000 (in μl) to use when preparing complexes should be 1:2 to 1:3 for most cell lines. To transfect $0.5 - 2 \times 10^5$ cells in a 24-well format, use 0.8-1 μg DNA and 2-3 μl of HY2000™. Optimizing transfection by varying DNA/HY2000 ratio is possible.
2. It is **CRITICAL** to transfect cells at high cell density. 90-95% confluence the time of transfection is recommended to obtain high efficiency and expression levels and to minimize decreased cell growth associated with high transfection activity. Lower cell densities are suitable with optimization of conditions. Take care to maintain a standard seeding protocol between experiments because transfection efficiency is dependent on culture confluence.
3. **DO NOT** add antibiotics to media during transfection as this will cause cell death.

For better results, you may choose to:

Use Opti-MEM I medium to dilute HY2000™ prior to complexing with DNA. Other media without serum (e.g. DMEM) may be used to dilute HY2000™ but transfection efficiency may be compromised.

Note: Some serum-free formulations can inhibit HY2000 mediated transfection, for example:

CD 293 · 293 SFM II · and VP-SFM etc.

Transfection Procedure for 24-Well Format

For adherent cells: One day before transfection · plate cells in growth medium (without antibiotics) so that they will be 90-95% confluent at the time of transfection ($0.5 - 2 \times 10^5$ cells/well for a 24-well plate).

For suspension cells: On the day of transfection just prior to preparing complexes · plate $4 - 8 \times 10^5$ cells/500 μl of growth medium (without antibiotics) in a 24-well plate.

1. For each transfection sample · prepare DNA-HY2000 complexes as follows:
 - Dilute DNA in 50 μl of Opti-MEM I Reduced Serum Medium without serum (or other medium without serum). Mix gently.
 - Mix HY2000 gently before use · then dilute the appropriate amount in 50 μl of Opti-MEM I Medium (or other medium without serum). Mix gently and incubate for 5 minutes at room temperature. Note: Combine the diluted HY2000 with the diluted DNA within 30 minutes. Longer incubation times may decrease activity. If DMEM is used as a diluent for the HY2000 · mix with the diluted DNA **within 5 minutes**.

- After the 5 minute incubation · combine the diluted DNA with the diluted HY2000™ (total volume is 100 µl). Mix gently and incubate for 20 minutes at room temperature to allow the DNA-HY2000 complexes to form. The solution may appear cloudy · but this will not inhibit the transfection. Note: DNA-HY2000 complexes are stable for at least 5 hours at room temperature.
3. Add the 100 µl of DNA-HY2000 complexes to each well. Mix gently by rocking the plate back and forth.
 4. Incubate the cells at 37°C in a CO₂ incubator for 24-48 hours until they are ready to assay for transgene expression. It is not necessary to remove the complexes or change the medium; however · growth medium may be replaced after 4-6 hours without loss of transfection activity.

For stable cell lines: Passage the cells at a 1:10 or higher dilution into fresh growth medium 24 hours after transfection. Add selective medium the following day.

For suspension cells: Add PMA and/or PHA (if desired) 4 hours after adding the DNA-HY2000 complexes to the cells. Tip: For Jurkat cells · adding PHA-L and PMA at final concentrations of 1 µg/ml and 50 ng/ml · respectively · enhances CMV promoter activity and gene expression. For K562 cells · adding PMA alone is sufficient to enhance promoter activity.

Scaling Up or Down Transfections

To transfect cells in different tissue culture formats · vary the amounts of HY2000 · DNA · cells · and medium used in proportion to the difference in surface area (see table below). With automated · high-throughput systems · larger complexing volumes are recommended for transfections in 96-well plates. Note: You may perform rapid 96-well plate transfections (plate cells and transfect simultaneously) by adding a suspension of cells directly to complexes prepared in the plate. Prepare complexes and add cells at twice the cell density as in the basic protocol in a 100 µl volume. Cells will adhere as usual in the presence of DNA-HY2000 complexes.

Culture Vessel	Surface Area per Well (cm ²)	Relative Surface Area (vs. 24-well)	Volume of Plating Medium	DNA (µg) and Dilution Volume (µl)	HY2000 (µl) and Dilution Volume (µl)
96-well	0.3	0.2	100 µl	0.2 µg in 25 µl	0.5 µl in 25 µl
24-well	2	1	500 µl	0.8 µg in 50 µl	2.0 µl in 50 µl
12-well	4	2	1 ml	1.6 µg in 100 µl	4.0 µl in 100 µl
35-mm	10	5	2 ml	4.0 µg in 250 µl	10 µl in 250 µl
6-well	10	5	2 ml	4.0 µg in 250 µl	10 µl in 250 µl
60-mm	20	10	5 ml	8.0 µg in 0.5 ml	20 µl in 0.5 ml
10-cm	60	30	15 ml	24 µg in 1.5 ml	60 µl in 1.5 ml

Note: Surface areas are determined from actual measurements of tissue culture vessels.

Optimizing Transfection

To obtain the highest transfection efficiency and low non-specific effects · optimize transfection conditions by varying DNA and HY2000 concentrations · and cell number. Make sure that cells are greater than 90% confluent and vary DNA (µg) : HY2000 (µl) ratios from 1:0.5 to 1:5.